

Root[®] with O3[®] Regional Oximetry

Available for Adult and Paediatric Applications



- > May help clinicians monitor cerebral oxygenation in situations in which peripheral pulse oximetry alone may not be fully indicative of the oxygen in the brain
- > Seamlessly integrates with the Root platform alongside SedLine[®] Brain Function Monitoring for a more complete brain monitoring solution

O3 Display

Δ base
Displays the difference between current rSO₂ and user-defined baseline

AUC
Area Under the Curve index quantifies the depth and duration of patient-stay below user-defined rSO₂ low alarm limit




rSO₂
Tissue oxygen saturation

rSO₂ Trend

Δ SpO₂
Displays the difference between SpO₂ (from the Radical-7®, if applicable) and rSO₂

O3 Monitoring

Root patient monitoring and connectivity hub offers plug and play monitoring with Masimo Open Connect™ (MOC-9™) modules.¹



Apply the appropriate O3 Sensors to the forehead:
 > Adult Adhesive Sensor (≥40 kg)
 > Paediatric Adhesive Sensor (<40 kg)



Connect the O3 Sensors to an O3 MOC-9 module (up to two sensors per module)



Connect the O3 MOC-9 module to one of three MOC-9 ports on Root

O3 MOC-9 Module Specifications

PHYSICAL CHARACTERISTICS	ENVIRONMENTAL
Length (including cable)..... 12.1 ft (3.7 m)	Operational Temperature..... 32 to 104°F (0 to 40°C)
Width..... 1.8 in (4.6 cm)	Storage Temperature..... -40 to 158°F (-40 to 70°C)
Thickness..... 0.6 in (1.5 cm)	Operating and Storage Humidity..... 10 to 95%, non-condensing
Weight..... 7.1 oz max (200 g max)	Altitude..... Up to 12,000 ft (3700 m)

O3 Series Sensor Specifications

Application Site..... Forehead	Wavelengths..... 4	ENVIRONMENTAL
Adult rSO₂ Sensor Accuracy (ARMS)² ≥40 kg	Absolute Regional Oxygen Saturation (rSO ₂)..... 4%	Operating Temperature at Ambient Humidity..... 41 to 104°F (5 to 40°C)
Trending Regional Oxygen Saturation (rSO ₂)..... 3%	Paediatric rSO₂ Sensor Accuracy (ARMS)² <40 kg	Storage Temperature at Ambient Humidity..... -40 to 140°F (-40 to 60°C)
Trending Regional Oxygen Saturation (rSO ₂)..... 3%		Storage Humidity..... 15 to 90%, 86 to 140°F (30 to 60°C)

¹In countries with regulatory approval and Root devices with the correct software version ² ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within ± ARMS of the reference measurements in a controlled study.

The O3 System with paediatric sensor is not licensed for sale in Canada.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

