

Masimo Trace™

Intuitive Patient Data Analysis and Reporting Tool



- > Create easy-to-read patient reports that include parameter trends, histograms, event annotations, and key statistics
- > Conveniently review reports for advanced Masimo measurements, including rainbow® parameters, brain monitoring, and capnography
- > Connect remotely to networked Masimo devices
- > Rapidly transfer up to 96 hours of device parameter data
- > Review and generate additional reports for past patient data
- > Facilitate clinician access to patient data at any time through network storage capability



Customisable and Easy-to-Read Reports

Generate a wide range of customisable reports that includes the following:



Additional Reporting Options

Respiratory Events Report
<ul style="list-style-type: none"> > Lists events according to a user-defined respiratory rate (RR) threshold and event duration > Lowest RR value is provided for each event
Events Report
<ul style="list-style-type: none"> > Lists user-added and parameter events such as low confidence, invalid measurements, and low and high alarms > Logs the duration and time of events
Time Interval Reports
<ul style="list-style-type: none"> > Graphically displays parameter measurements over 24, 12, 8, 2 and 1 hour and 30, 10 and 1 minute time intervals

Easy Installation

1. Plug the Trace USB drive into a USB port on the Trace host computer.
2. Locate the Trace installation (.exe) file.
3. Double click the Trace Setup.exe file to start the Trace installation process.
4. Read and accept the Masimo end user license agreement.
5. If a Desktop Shortcut to Trace is desired, check the box next to "Create a desktop shortcut". Click Next to continue.
6. Click Next to perform the installation.
7. By default, the check box to "Launch Masimo Trace" is selected (de-select if desired). Once Trace is installed, click Finish. Trace will open.

Technical Information

PC SOFTWARE REQUIREMENTS	REQUIRED CABLES
Trace is compatible with PCs running Windows 7 or higher	Root and Rad-97 USB-Null Modem Cable
SUPPORTED DEVICES	Radical-7 (docked to RDS) USB-to-RS-232 Serial Adapter Cable
Radical-7® (V1.4.5.1 or higher) with RDS (V5.1.3.2 or higher)	INCLUDED IN TRACE REPORTS
Root® (V16.2.1 or higher) with Radical-7 (V1.4.5.1 or higher) and connected MOC-9™ modules	SET®
Root (V16.2.1 or higher) with Radius-7® and connected MOC-9 modules	rainbow®
Rad-97™ product family (V1.0.3.5 or higher)	Respiration Rate
CONNECTION TYPES FOR DATA DOWNLOAD	Capnography
USB or USB to RS-232 Serial Cable (depending on Masimo device)	SedLine® Brain Function Monitoring
Network, using the Masimo device's IP address ³	O3® Regional Oximetry
ORDERING INFORMATION	Noninvasive Blood Pressure and Temperature
Masimo Trace USB Installer. PN 9853	FILES GENERATED
LICENSE	CSV or PDF
Valid for one year	

Available in English only

¹ Available when SpO2 parameter is selected for inclusion in reports.

² Available when Acoustic Respiration Rate (RRa®) is selected for inclusion in reports.

³ Trace does not support network connections for Rad-97.