

# Kite™

Expanded Visibility of Patient Data




Kite provides a supplemental display of data from a Masimo device


- > Display patient data from Masimo monitors on compatible smart devices
- > Configure the display differently than that of the connected Masimo device
- > Integrates into existing hospital infrastructures where a supplemental display may be beneficial, such as the operating room or cardiac theater

## Configurable Display Options


By customising what is displayed, Kite allows clinicians to focus on the most pertinent data for each stage of a patient's care, empowering clinicians to make more informed decisions



Configure Kite's supplemental display to feature the most applicable data for the current patient or case



Display data from modules connected through MOC-9™ ports to Root®, such as SedLine® brain function monitoring, O3™ regional oximetry, or an ISA™ EtCO2 module



Displays visual alarms from the patient monitor, providing quick notification of changes in a patient's physiologic status

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

**Masimo U.S.**  
Tel: 1 877 4 Masimo  
info-america@masimo.com

**Masimo International**  
Tel: +41 32 720 1111  
info-international@masimo.com

