





NomoLine® Catalog




Moisture-wicking sampling lines for intubated and non-intubated patients in low- and high-humidity applications



- > Patented moisture-wicking technology allows water evaporation from the sampling line
- > Hydrophobic bacterial filter prevents contamination of the gas analyser
- > Uniquely designed sampling line allows for usage during low-flow and small sample-size applications
- > Ergonomically designed cannulas improve patient comfort
- > Star Lumen O2 delivery tubing allows for continuous, accurate O2 supply even if kinked








NomoLine Sampling Line Options






Division	Product Description	Part Number	Length	Image
Low Humidity or Short-term Airway Adapter Set	NomoLine LH Adult/Paediatric Airway Adapter Set, 25/box	3814	2m	
	NomoLine LH Adult/Paediatric Airway Adapter Set, 4m, 25/box	4620	4m	
High Humidity or Long-term Airway Adapter Set	NomoLine HH Adult/Paediatric Airway Adapter Set, 25/box	3827	2m	
	NomoLine HH Adult/Paediatric Airway Adapter Set, 25/box	3828	3m	
	NomoLine HH Infant/Neonatal, Airway Adapter Set, 25/box	4367	2m	
Low Humidity or Short-term Cannula with O2 Tubing	NomoLine LH Adult Nasal CO2 Cannula with O2, 25/box	3820	2m	
	NomoLine LH Adult Nasal CO2 Cannula with O2, 4m, 25/box	4624	4m	
	NomoLine LH Paediatric Nasal CO2 Cannula with O2, 25/box	3821	2m	
	NomoLine LH Adult Nasal/Oral CO2 Cannula with O2, 25/box	3824	2m	
	NomoLine LH Adult Nasal/Oral CO2 Cannula with O2, 4m, 25/box	4628	4m	
	NomoLine LH Paediatric Nasal/Oral CO2 Cannula with O2, 25/box	3825	2m	
	NomoLine LH Paediatric Nasal/Oral CO2 Cannula with O2, 4m, 25/box	4629	4m	
	NomoLine HH Adult Nasal CO2 Cannula with O2, 25/box	3833	2m	
NomoLine HH Paediatric Nasal CO2 Cannula with O2, 25/box	3834	2m		
NomoLine HH Adult Nasal/Oral CO2 Cannula with O2, 25/box	3837	2m		
NomoLine HH Paediatric Nasal/Oral CO2 Cannula with O2, 25/box	3838	2m		

Division	Product Description	Part Number	Length	Image	
Low Humidity or Short-term Cannula without O2 tubing	NomoLine LH Adult Nasal CO2 Cannula, 25/box	3817	2m		
	NomoLine LH Paediatric Nasal CO2 Cannula, 25/box	3818	2m		
	NomoLine LH Infant Nasal CO2 Cannula, 25/box	3819	2m		
	High Humidity or Long-term Cannula without O2 tubing	NomoLine LH Adult Nasal/Oral CO2 Cannula, 25/box	3822	2m	
		NomoLine LH Adult Nasal/Oral CO2 Cannula, 4m, 25/box	4626	4m	
		NomoLine LH Paediatric Nasal/Oral CO2 Cannula, 25/box	3823	2m	
		NomoLine LH Adult Single Nasal Prong CO2 Cannula, 25/box	3826	2m	
		NomoLine HH Adult Nasal CO2 Cannula, 25/box	3830	2m	
NomoLine HH Infant Nasal CO2 Cannula, 25/box	3832	2m			
NomoLine HH Adult Single Nasal Prong CO2 Cannula, 25/box	3839	2m			

LH = Low Humidity HH = High Humidity

NomoLine-O™ Sampling Line Options

Division	Product Description	Part Number	Length	Image
Low Humidity or Short-term Airway Adapter Set	NomoLine-O LH Adult/Paediatric Airway Adapter Set, 2m, 25/Box	4454	2m	
	NomoLine-O LH Adult/Paediatric Airway Adapter Set, 4m, 25/Box	4630	4m	
High Humidity or Long-term Airway Adapter Set	NomoLine-O HH Adult/Paediatric Airway Adapter Set, 2m, 25/Box	4449	2m	
	NomoLine-O HH Infant/Neonatal Airway Adapter Set, 2m, 25/Box	4450	2m	
Low Humidity or Short-term Cannula with O2 Tubing	NomoLine-O LH Adult Nasal CO2 Cannula with O2, 2m, 25/Box	4456	2m	
	NomoLine-O LH Adult Nasal CO2 Cannula with O2, 4m, 25/Box	4634	4m	
	NomoLine-O LH Paediatric Nasal CO2 Cannula with O2, 2m, 25/Box	4457	2m	
	NomoLine-O LH Adult Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4459	2m	
	NomoLine-O LH Adult Nasal/Oral CO2 Cannula with O2, 4m, 25/Box	4638	4m	
	NomoLine-O LH Paediatric Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4460	2m	
	NomoLine-O LH Paediatric Nasal/Oral CO2 Cannula with O2, 4m, 25/Box	4639	4m	
	NomoLine-O LH Paediatric Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4460	2m	
High Humidity or Long-term Cannula with O2 tubing	NomoLine-O HH Adult Nasal CO2 Cannula with O2, 2m, 25/Box	4451	2m	
	NomoLine-O HH Paediatric Nasal CO2 Cannula with O2, 2m, 25/Box	4645	2m	
	NomoLine-O HH Adult Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4453	2m	
	NomoLine-O HH Paediatric Nasal/Oral CO2 Cannula with O2, 2m, 25/Box	4647	2m	

Division	Product Description	Part Number	Length	Image	
Low Humidity or Short-term Cannula without O2 Tubing	NomoLine-O LH Adult Nasal CO2 Cannula, 2m, 25/box	4455	2m		
	NomoLine-O LH Paediatric Nasal CO2 Cannula, 2m, 25/box	4641	2m		
	NomoLine-O LH Infant Nasal CO2 Cannula, 2m, 25/box	4662	2m		
	High Humidity or Long-term Cannula without O2 tubing	NomoLine-O LH Adult Nasal/Oral CO2 Cannula, 2m, 25/Box	4458	2m	
		NomoLine-O LH Adult Nasal/Oral CO2 Cannula, 4m, 25/Box	4636	4m	
		NomoLine-O LH Paediatric Nasal/Oral CO2 Cannula, 2m, 25/box	4642	2m	
		NomoLine-O LH Adult Single Nasal Prong CO2 Cannula, 2m, 25/box	4663	2m	
High Humidity or Long-term Cannula without O2 tubing	NomoLine-O HH Adult Nasal CO2 Cannula, 2m, 25/box	4643	2m		
	NomoLine-O HH Infant Nasal CO2 Cannula, 2m, 25/box	4664	2m		
	NomoLine-O HH Adult Single Nasal Prong CO2 Cannula, 2m, 25/box	4665	2m		

LH = Low Humidity HH = High Humidity

NomoLine with Luer Sampling Line Options

Division	Product Description	Part Number	Length	Image
Airway Adapter Set	NomoLine Adult/Paediatric Airway Adapter Set with Luer, 2m, 25/box	3801	2m	
Cannula with O2 Tubing	NomoLine Adult Nasal CO2 Cannula with O2 and Luer, 2m, 25/box	3807	2m	
	NomoLine Paediatric Nasal CO2 Cannula with O2 and Luer, 2m, 25/box	3808	2m	
	NomoLine Adult Nasal/Oral CO2 Cannula with O2 and Luer, 2m, 25/box	3811	2m	
	NomoLine Paediatric Nasal/Oral CO2 Cannula with O2 and Luer, 2m, 25/box	3812	2m	
Cannula without O2 Tubing	NomoLine Adult Nasal/Oral CO2 Cannula with Luer, 2m, 25/box	3809	2m	
	NomoLine Adult Single Nasal Prong CO2 Cannula with Luer, 2m, 25/box	3813	2m	

Sampling Line Connector Versions



HH Sampling Lines - ISA™

Long-term usage

For use with Masimo ISA and Rad-97® devices



LH Sampling Lines - ISA

Short-term usage

For use with Masimo ISA and Rad-97 devices



HH Sampling Lines - Medtronic® (Oridion)

Long-term usage

For use with Medtronic Microstream™ devices and modules



LH Sampling Lines - Medtronic (Oridion)

Short-term usage

For use with Medtronic Microstream devices and modules



Luer Sampling Lines - Water Trap

For use with luer-based devices using water traps

Not all part numbers are available for sale in Canada.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions

The Medtronic® trademarks are property of Medtronic and are used for identification purposes only. Use of the name, trademark, or brand does not imply any affiliation or endorsement.

Masimo Sweden AB
Tel: +46 8 544 98 150
info-sweden@masimo.com

Masimo U.S.
tel 1-877-4-Masimo
info-america@masimo.com

Masimo International
tel +41-32-720-1111
info-international@masimo.com

