

# UniView™

Intelligent Display for Next Generation Data Aggregation and Alarm Visualisation




- > UniView aggregates data and alarms from multiple Masimo and third-party devices – such as patient monitors, ventilators, anaesthesia machines, IV pumps and others connected through Masimo systems – on a supplemental display
- > Integrated real-time data visualisation reduces cognitive overload and promotes data sharing among multiple clinicians, enabling them to spot trends and coordinate care
- > Visual alarm indicators, relayed from connected devices, help care teams recognise patient distress and target areas for clinical focus
- > Tailored use-case-specific screen layouts optimise the presentation of advanced and integrated parameters, trend data, and waveforms in critical care areas
- > Adaptable layout automatically reconfigures based on connected devices

## Tailored Screen Layouts


Use-case-specific screen layouts allow clinicians to view the most pertinent data applicable to a clinical care team

### Overview Layout




Adaptable layout automatically reconfigures based on connected devices to optimise the view of parameter data, including waveforms and alarms, for an overview of patient status

### Haemodynamics Layout



View trend data for noninvasive haemoglobin (SpHb®), pleth variability index (PVI®), and pulse rate (PR) to aid in visualising a patient's status over time

### Oxygenation Layout



View ventilator waveforms alongside noninvasive trended haemoglobin (SpHb) and oxygen saturation (SpO2) to monitor a patient's oxygenation status

### Sedation Layout



View high-fidelity EEG waveforms, patient state index (PSI), and anaesthesia machine data to monitor a patient's depth of anaesthesia

## Ordering Information

	Part Number
UniView Mini-PC .....	9856
<b>Required Masimo System Components</b>	
Patient SafetyNet™* .....	9095
Iris Gateway™ .....	9712

\*The use of the trademark SafetyNet is under license from University HealthSystem Consortium.

ORi™ has obtained CE Marking. Not available in the U.S.  
ORi is not licensed for sale in Canada.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

**Masimo U.S.**  
Tel: 1 877 4 Masimo  
info-america@masimo.com

**Masimo International**  
Tel: +41 32 720 1111  
info-international@masimo.com

